




CT455 rabbit pAb

Catalog No	YP-Ab-12307
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CT45A5 CT45-5
Immunogen	Synthesized peptide derived from human CT455 AA range: 135-185
Specificity	This antibody detects endogenous levels of CT455 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CT455;CT45A5 CT45-5
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	tissue specificity:Testis specific. Expressed in cancer cell lines.,
Background	This gene represents one of a cluster of several similar genes located on the q arm of chromosome X. The genes in this cluster encode members of the cancer/testis (CT) family of antigens, and are distinct from other CT antigens. These antigens are thought to be novel therapeutic targets for human cancers. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4° over night