




LARP4 rabbit pAb


Catalog No	YP-Ab-12309
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;ELISA;IHC
Gene Name	LARP4 PP13296
Immunogen	Synthesized peptide derived from human LARP4 AA range: 209-259
Specificity	This antibody detects endogenous levels of LARP4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	La related protein 4;LARP4;PP13296;LARP4 PP13296
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, Stress granule . Cytoplasm, cytosol . Localized throughout the cytosol. Partially localized in stress granules in response to arsenite treatment. .
Tissue Specificity	
Function	sequence caution:Translation N-terminally extended.,similarity:Contains 1 HTH La-type RNA-binding domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).