



PIWL4 rabbit pAb

Catalog No	YP-Ab-12315
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; IHC
Gene Name	PIWIL4 HIWI2 PIWI
Immunogen	Synthesized peptide derived from human PIWL4 AA range: 108-158
Specificity	This antibody detects endogenous levels of PIWL4 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DKFZp686P1248; FLJ36156; HILI 2; HILI2; HIWI 2; HIWI2; Miwi 2 protein; Miwi2; PIWI; Piwi like 2; Piwi like 4 (Drosophila); Piwi like 4; Piwi like protein 4; PIWI like protein; Piwi like RNA mediated gene silencing 4; Piwi-like protein 4; PIWIL 4; Piwi4; PIWL4_HUMAN; PIWL4; PIWIL4 HIWI2 PIWI
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Cytoplasm . Probable component of the meiotic nuage, also named P granule, a germ-cell-specific organelle required to repress transposon activity during meiosis. PIWIL2/MILI is required for nuclear localization (By similarity).
Tissue Specificity	Ubiquitously expressed (PubMed:25038252, PubMed:17544373, PubMed:28025795, PubMed:28711973, PubMed:22483988). Detected in retina, retinal pigment epithelia cells (RPE) (at protein level) (PubMed:28025795).
Function	function: Plays a central role during spermatogenesis by repressing transposable elements and prevent their mobilization, which is essential for the germline integrity. Acts via the piRNA metabolic process, which mediates the repression of transposable elements during meiosis by forming complexes composed of piRNAs and Piwi proteins and govern the methylation and subsequent repression of transposons. Directly binds piRNAs, a class of 24 to 30 nucleotide



RNAs that are generated by a Dicer-independent mechanism and are primarily derived from transposons and other repeated sequence elements. Associates with secondary piRNAs antisense and PIWIL2/MILI is required for such association. The piRNA process acts upstream of known mediators of DNA methylation. Participates to a piRNA amplification loop. Besides their function in transposable elements repression, piRNAs are probably involved in other

Background

PIWIL4 belongs to the Argonaute family of proteins, which function in development and maintenance of germline stem cells (Sasaki et al., 2003 [PubMed 12906857]).[supplied by OMIM, Mar 2008],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human cervical carcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).