




KANK3 rabbit pAb

Catalog No	YP-Ab-12316
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	KANK3 ANKRD47
Immunogen	Synthesized peptide derived from human KANK3 AA range: 135-185
Specificity	This antibody detects endogenous levels of KANK3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KANK3;KANK3 ANKRD47
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Strongly expressed in breast, liver, lung, skeletal muscle and kidney.
Function	similarity:Contains 5 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night