




TECT3 rabbit pAb

Catalog No	YP-Ab-12320
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TCTN3 C10orf61 TECT3 PSEC0041 UNQ1881/PRO4324
Immunogen	Synthesized peptide derived from human TECT3 AA range: 497-547
Specificity	This antibody detects endogenous levels of TECT3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TECT3;TCTN 3;PSEC0041;OFD4;JBTS18;TCTN3 C10orf61 TECT3 PSEC0041 UNQ1881/PRO4324
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	function:May be involved in apoptosis regulation.,similarity:Belongs to the tectonic family.,
Background	This gene encodes a member of the tectonic gene family which functions in Hedgehog signal transduction and development of the neural tube. Mutations in this gene have been associated with Orofaciodigital Syndrome IV and Joubert Syndrom 18. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2012],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from CACO2 cells, primary antibody was diluted at 1:1000, 4° over night