




TEPP rabbit pAb

Catalog No	YP-Ab-12327
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	TEPP
Immunogen	Synthesized peptide derived from human TEPP AA range: 24-74
Specificity	This antibody detects endogenous levels of TEPP at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TEPP
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed in testis, prostate and placenta.
Function	tissue specificity:Expressed in testis, prostate and placenta.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night