




AGAP9 rabbit pAb

Catalog No	YP-Ab-12342
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	AGAP9 CTGLF6
Immunogen	Synthesized peptide derived from human AGAP9 AA range: 4-54
Specificity	This antibody detects endogenous levels of AGAP9 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AGAP9;AGAP9 CTGLF6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., function:Putative GTPase-activating protein .,miscellaneous:Encoded by one of the numerous copies of centaurin gamma-like genes clustered in the q11 region of chromosome 10.,similarity:Belongs to the centaurin gamma-like family., similarity:Contains 1 ANK repeat.,similarity:Contains 1 Arf-GAP domain., similarity:Contains 1 PH domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night