




UNKL rabbit pAb

Catalog No	YP-Ab-12346
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	UNKL C16orf28 ZC3H5L ZC3HDC5L
Immunogen	Synthesized peptide derived from human UNKL AA range: 296-346
Specificity	This antibody detects endogenous levels of UNKL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C16orf28;EC:2.3.2.-;KIAA0734;Putative E3 ubiquitin-protein ligase UNKL;RING finger protein unkempt-like;UNKL;UNKL C16orf28 ZC3H5L ZC3HDC5L
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 4]: Cytoplasm. Nucleus. Isoform 4 is primarily localized in the cytoplasm but has the ability to shuttle between the nucleus and the cytoplasm.
Tissue Specificity	
Function	similarity:Belongs to the unkempt family.,similarity:Contains 1 RING-type zinc finger.,
Background	This gene encodes a RING finger protein that may function in Rac signaling. It can bind to Brg/Brm-associated factor 60b and can promote its ubiquitination. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2013],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night