




T194A rabbit pAb

Catalog No	YP-Ab-12352
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM194A KIAA0286 TMEM194
Immunogen	Synthesized peptide derived from human T194A AA range: 164-214
Specificity	This antibody detects endogenous levels of T194A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0286;TMEM194;TMEM194A;transmembrane protein 194A;T194A; TMEM194A KIAA0286 TMEM194
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus inner membrane ; Multi-pass membrane protein ; Nucleoplasmic side . Nucleus envelope . Colocalizes with lamins and RAN-GTP at the nuclear envelope. .
Tissue Specificity	
Function	similarity:Belongs to the TMEM194 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night