




NET5 rabbit pAb

Catalog No	YP-Ab-12365
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NTN5
Immunogen	Synthesized peptide derived from human NET5 AA range: 91-141
Specificity	This antibody detects endogenous levels of NET5 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NET 5;NET5;PP1057;Tetraspan NET 5;tetraspanin 9;Tspan 9;TSPAN9;NTN5
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity:Contains 1 NTR domain.,similarity:Contains 3 laminin EGF-like domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night