




F90A9 rabbit pAb

Catalog No	YP-Ab-12366
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM90A9P FAM90A9
Immunogen	Synthesized peptide derived from human F90A9 AA range: 70-120
Specificity	This antibody detects endogenous levels of F90A9 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F90A9;FAM90A9P FAM90A9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	caution:Could be the product of a pseudogene.,similarity:Belongs to the FAM90 family.,
Background	FAM90A9 belongs to subfamily II of the primate-specific FAM90A gene family, which originated from multiple duplications and rearrangements (Bosch et al., 2007 [PubMed 17684299]). For background information on the FAM90A gene family, as well as information on the evolution of FAM90A genes, see FAM90A1 (MIM 613041).[supplied by OMIM, Oct 2009],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night