




ZIK1 rabbit pAb

Catalog No	YP-Ab-12379
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZIK1 ZNF762
Immunogen	Synthesized peptide derived from human ZIK1 AA range: 171-221
Specificity	This antibody detects endogenous levels of ZIK1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZIK1;Zinc finger protein 762;ZNF762;ZIK1 ZNF762
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed at high levels in gastric glands, and at low levels in colon and small intestine. Silenced through promoter methylation in gastric glands with intestinal metaplasia.
Function	function:May be a transcriptional repressor.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 9 C2H2-type zinc fingers.,subunit:Interacts with HNRPK.,tissue specificity:Expressed at high levels in gastric glands, and at low levels in colon and small intestine. Silenced through promoter methylation in gastric glands with intestinal metaplasia.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night