




F209A rabbit pAb

Catalog No	YP-Ab-12383
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM209A C20orf106
Immunogen	Synthesized peptide derived from human F209A AA range: 65-115
Specificity	This antibody detects endogenous levels of F209A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F209A;FAM209A C20orf106
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus inner membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	sequence caution:Contaminating sequence. Potential poly-A sequence.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night