




ATD3C rabbit pAb

Catalog No	YP-Ab-12384
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ATAD3C
Immunogen	Synthesized peptide derived from human ATD3C AA range: 354-404
Specificity	This antibody detects endogenous levels of ATD3C at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ATD3C;ATAD3C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,sequence caution:Wrong choice of frame.,similarity:Belongs to the AAA ATPase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night