




SKDA1 rabbit pAb

Catalog No	YP-Ab-12385
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SKIDA1 C10orf140 DLN1
Immunogen	Synthesized peptide derived from human SKDA1 AA range: 714-764
Specificity	This antibody detects endogenous levels of SKDA1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SKDA1;SKIDA1 C10orf140 DLN1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	polymorphism:The poly-Ala region of C10orf140 is highly polymorphic and the number of Ala can vary.,similarity:Belongs to the DACH/dachshund family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night