




BRD9 rabbit pAb

Catalog No	YP-Ab-12390
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	BRD9 UNQ3040/PRO9856
Immunogen	Synthesized peptide derived from human BRD9 AA range: 310-360
Specificity	This antibody detects endogenous levels of BRD9 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BRD 9;Bromodomain-containing protein 9;LAVS3040;Rhabdomyosarcoma antigen MU-RMS-40.8;UNQ3040/PRO9856;BRD9;BRD9 UNQ3040/PRO9856
Observed Band	66kD
Calculated Molecular Weight	66kD
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 bromo domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night