



KITH_EBV rabbit pAb

Catalog No	YP-Ab-12409
Isotype	IgG
Reactivity	0
Applications	WB; ELISA
Gene Name	KITH
Immunogen	Synthesized peptide derived from KITH_EBV
Specificity	This antibody detects endogenous levels of KITH_EBV
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KITH_EBV;KITH
Observed Band	
Calculated Molecular Weight	67kDa;
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

