



KR2_HHV11 rabbit pAb

Catalog No	YP-Ab-12410
Isotype	IgG
Reactivity	0
Applications	WB; ELISA
Gene Name	KR2
Immunogen	Synthesized peptide derived from KR2_HHV11
Specificity	This antibody detects endogenous levels of KR2_HHV11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KR2_HHV11;KR2
Observed Band	
Calculated Molecular Weight	57 kDa;
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

