



ALS2CR16 rabbit pAb

Catalog No	YP-Ab-12415
Isotype	IgG
Reactivity	0
Applications	WB; ELISA
Gene Name	ALS2CR16
Immunogen	Synthesized peptide derived from ALS2CR16
Specificity	This antibody detects endogenous levels of ALS2CR16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALS2CR16;ALS2CR17;neurobeachin like 1;Neurobeachin like protein 1
Observed Band	280 kDa
Calculated Molecular Weight	200 aa, 307 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

