



GNG5 rabbit pAb

Catalog No	YP-Ab-12419
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB; ELISA
Gene Name	GNG5
Immunogen	Synthesized peptide derived from human GNG5
Specificity	This antibody detects endogenous levels of Human GNG5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ92393;GBG5_HUMAN;Gng5;Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-5;OTTHUMP11474;OTTHUMP11565
Observed Band	
Calculated Molecular Weight	8 kDa; 7kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

