



GFI1B rabbit pAb

Catalog No	YP-Ab-12423
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB; ELISA
Gene Name	GFI1B
Immunogen	Synthesized peptide derived from human GFI1B
Specificity	This antibody detects endogenous levels of Human GFI1B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GFI 1B;gfi1b;GFI1B protein;GFI1B_HUMAN;Growth factor independent 1B protein;Growth factor independent 1B transcription repressor;Growth factor independent protein 1B;OTTHUMP22443;OTTHUMP22444;OTTHUMP235527; Potential regulator of CDKN1A;Potential regulator of CDKN1A translocated in CML;Translocated in CML;Zinc finger protein Gfi-1b
Observed Band	
Calculated Molecular Weight	43 kDa; 37kD(Calculated).
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

