



DUET rabbit pAb

Catalog No	YP-Ab-12425
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB; ELISA
Gene Name	DUET
Immunogen	Synthesized peptide derived from human DUET
Specificity	This antibody detects endogenous levels of Human DUET
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DUET;duo;HAPIP;Kalirin;RhoGEF kinase;KALRN;Protein Duo;TRAD
Observed Band	200 kDa
Calculated Molecular Weight	340 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

