



IL18BP rabbit pAb

Catalog No	YP-Ab-12426
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB; ELISA
Gene Name	IL18BP
Immunogen	Synthesized peptide derived from human IL18BP AA range: 120-200
Specificity	This antibody detects endogenous levels of Human IL18BP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IL18BP;IL 18BP;IL18BPα;interleukin 18 binding protein;Interleukin-18-binding protein
Observed Band	21 kDa
Calculated Molecular Weight	194 aa, 21 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

