



AGRP rabbit pAb

Catalog No	YP-Ab-12441
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB; ELISA
Gene Name	AGRP AGRT ART
Immunogen	Synthesized peptide derived from human AGRP AA range: 21-70
Specificity	This antibody detects endogenous levels of Human AGRP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Agouti related protein;AGRP;AGRT;ART;ASIP2;KNOTTIN;AGRP AGRT ART; Agouti-related protein
Observed Band	
Calculated Molecular Weight	132 aa, 14 kDa
Cell Pathway	Secreted . Golgi apparatus lumen .
Tissue Specificity	Expressed primarily in the adrenal gland, subthalamic nucleus, and hypothalamus, with a lower level of expression occurring in testis, lung, and kidney.
Function	cell surface receptor linked signal transduction, G-protein coupled receptor protein signaling pathway, neuropeptide signaling pathway, intracellular signaling cascade, behavior, feeding behavior, adult feeding behavior, response to endogenous stimulus, response to hormone stimulus, hormone-mediated signaling, response to organic substance,adult behavior, cellular response to hormone stimulus, eating behavior,
Background	disease:Defects in AGRP may be a cause of autosomal dominant obesity [MIM:601665].,domain:The presence of a 'disulfide through disulfide knot' structurally defines this protein as a knottin.,function:Plays a role in weight

homeostasis. May play a role in the regulation of melanocortin receptors within the hypothalamus and adrenal gland, and therefore in the central control of feeding.,similarity:Contains 1 agouti domain.,tissue specificity:Expressed primarily in the adrenal gland, subthalamic nucleus, and hypothalamus, with a lower level of expression occurring in testis, lung, and kidney.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images