



IL25(IL17E) rabbit pAb

Catalog No	YP-Ab-12467
Isotype	IgG
Reactivity	Mouse
Applications	WB; ELISA
Gene Name	IL25,IL-17E
Immunogen	Synthesized peptide derived from mouse IL25(IL17E)
Specificity	This antibody detects endogenous levels of Mouse IL17E
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IL17E;IL25;IL25(IL17E);IL-17E;Interleukin-25;Interleukin-17E
Observed Band	20kDa
Calculated Molecular Weight	20kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

