



IL27B rabbit pAb

Catalog No	YP-Ab-12468
Isotype	IgG
Reactivity	Mouse
Applications	WB; ELISA
Gene Name	IL27B
Immunogen	Synthesized peptide derived from mouse IL27B
Specificity	This antibody detects endogenous levels of Mouse IL27B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IL27B;IL35B;IL-27B;EBI3;Interleukin-27 subunit beta (IL-27 subunit beta;Epstein-Barr virus-induced gene 3 protein;EBV-induced gene 3 protein)
Observed Band	31kDa
Calculated Molecular Weight	25kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

