



IL5 rabbit pAb

Catalog No	YP-Ab-12469
Isotype	IgG
Reactivity	Mouse
Applications	WB; ELISA
Gene Name	IL5
Immunogen	Synthesized peptide derived from C-terminal mouse IL5
Specificity	This antibody detects endogenous levels of Mouse IL5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EDF; IL 5; IL5; IL-5; Interleukin 5; T cell replacing factor; TRF; Interleukin-5 (IL-5; B-cell differentiation factor I; Eosinophil differentiation factor; T-cell replacing factor; TRF)
Observed Band	
Calculated Molecular Weight	134 aa, 15 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

