



FGF-7 rabbit pAb

Catalog No	YP-Ab-12487
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	FGF7 KGF
Immunogen	Synthesized peptide derived from human FGF-7
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat FGF-7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FGF7;FGF-7;FGF 7;Fibroblast growth factor 7;HBGF-7;FGF7 KGF;Fibroblast growth factor 7 (FGF-7;Heparin-binding growth factor 7;Keratinocyte growth factor)
Observed Band	23 kDa
Calculated Molecular Weight	22-28 kDa
Cell Pathway	Secreted.
Tissue Specificity	Epithelial cell.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

