




POU4F2 rabbit pAb

Catalog No	YP-Ab-12497
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA;IHC
Gene Name	POU4F2 BRN3B
Immunogen	Synthesized peptide derived from human POU4F2
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat POU4F2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Brn-3B;BRN3B;Brain-specific homeobox/POU domain protein 3B;Brain-3B;POU4F2;POU4F2 BRN3B;POU domain;class 4;transcription factor 2 (Brain-specific homeobox/POU domain protein 3B;Brn-3B)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Nucleus speckle . Cytoplasm .
Tissue Specificity	Expressed in the brain (PubMed:7691107). Expressed in the ganglion cell layer of the retina (PubMed:7691107).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).