



A0A3Q3A7D8 rabbit pAb

Catalog No	YP-Ab-12498
Isotype	IgG
Reactivity	Chicken
Applications	WB; ELISA
Gene Name	A0A3Q3A7D8
Immunogen	Synthesized peptide derived from Chicken A0A3Q3A7D8
Specificity	This antibody detects endogenous levels of Chicken A0A3Q3A7D8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A0A3Q3A7D8
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

