



NAGK rabbit pAb

Catalog No	YP-Ab-12502
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	NAGK
Immunogen	Synthesized peptide derived from human NAGK
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat NAGK
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSA242910;GNK;GlcNAc kinase;EC:2.7.1.60;EC:2.7.1.59;NAGK;N-acetyl-D-glucosamine kinase (N-acetylglucosamine kinase;EC 2.7.1.59;GlcNAc kinase)
Observed Band	37 kDa
Calculated Molecular Weight	37 kDa
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

