



IDNK rabbit pAb

Catalog No	YP-Ab-12504
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	IDNK C9orf103
Immunogen	Synthesized peptide derived from human IDNK
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat IDNK
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C9orf103;EC:2.7.1.12;IDNK;IDNK C9orf103;Probable gluconokinase (EC 2.7.1.12;Gluconate kinase)
Observed Band	21 kDa
Calculated Molecular Weight	230 aa, 25 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

