



IL-3 rabbit pAb

Catalog No	YP-Ab-12513
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	IL3
Immunogen	Synthesized peptide derived from human IL-3
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat IL-3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Hematopoietic growth factor;IL 3;IL3;IL-3;Interleukin 3;Mast cell growth factor; MCGF;MULTI CSF;P cell stimulating factor;Interleukin-3 (IL-3;Multipotential colony-stimulating factor;P-cell-stimulating factor)
Observed Band	20 kDa
Calculated Molecular Weight	152 aa, 17 kDa
Cell Pathway	Secreted.
Tissue Specificity	Activated T-cells, mast cells, natural killer cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

