



TET3 rabbit pAb

Catalog No	YP-Ab-12534
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	TET3 KIAA0401
Immunogen	Synthesized peptide derived from human TET3
Specificity	This antibody detects endogenous levels of TET3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.16% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	tet oncogene family member 3;Methylcytosine dioxygenase TET3;KIAA0401; EC:1.14.11.80;TET3;TET3 KIAA0401;Methylcytosine dioxygenase TET3 (EC 1.14.11.n2)
Observed Band	180 kDa
Calculated Molecular Weight	1660 aa, 179 kDa
Cell Pathway	Nucleus . Cytoplasm . Chromosome . At the zygotic stage, localizes in the male pronucleus, while it localizes to the cytoplasm at other preimplantation stages. Binds to the promoter of target genes, close to the transcription start site. .
Tissue Specificity	Expressed in colon, muscle, adrenal gland and peripheral blood lymphocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

