



PSMC1 rabbit pAb

Catalog No	YP-Ab-12575
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	PSMC1
Immunogen	Synthesized peptide derived from human PSMC1
Specificity	This antibody detects endogenous levels of PSMC1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.66% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	26S proteasome AAA-ATPase subunit RPT2;26S proteasome regulatory subunit 4;P26S4;PSMC1;26S protease regulatory subunit 4 (P26s4) (26S proteasome AAA-ATPase subunit RPT2) (Proteasome 26S subunit ATPase 1)
Observed Band	52-57 kDa
Calculated Molecular Weight	49 kDa
Cell Pathway	Cytoplasm. Nucleus. Membrane ; Lipid-anchor .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

