



tet1 rabbit pAb

Catalog No	YP-Ab-12580
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	TET1
Immunogen	Synthesized peptide derived from human tet1
Specificity	This antibody detects endogenous levels of tet1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.72% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LCX;CXXC6;bA119F7.1;TET1
Observed Band	310kDa
Calculated Molecular Weight	235kDa
Cell Pathway	Nucleus . Chromosome .
Tissue Specificity	Expressed in fetal heart, lung and brain, and in adult skeletal muscle, thymus and ovary. Not detected in adult heart, lung or brain. Up-regulated in glioblastoma cells (at protein level) (PubMed:25284789).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

