




KLHL11 Monoclonal Antibody


Catalog No	YP-Ab-12605
Isotype	IgG
Reactivity	Human
Applications	WB;IF;FCM;ELISA
Gene Name	KLHL11
Immunogen	Purified recombinant fragment of human KLHL11 expressed in E. Coli.
Specificity	KLHL11 Monoclonal Antibody detects endogenous levels of KLHL11 protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Kelch-like protein 11; KLHL11
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm, SCF ubiquitin ligase complex,
Tissue Specificity	Testis,
Function	similarity: Contains 1 BACK (BTB/Kelch associated) domain., similarity: Contains 1 BTB (POZ) domain., similarity: Contains 5 Kelch repeats.,
Background	similarity: Contains 1 BACK (BTB/Kelch associated) domain., similarity: Contains 1 BTB (POZ) domain., similarity: Contains 5 Kelch repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Western Blot analysis using KLHL11 Monoclonal Antibody against HeLa (1) and MCF-7 (2) cell lysate.



Immunofluorescence analysis of 3T3-L1 cells using KLHL11 Monoclonal Antibody (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of K562 cells using KLHL11 Monoclonal Antibody (green) and negative control (purple).