




KLHL25 Monoclonal Antibody

Catalog No	YP-Ab-12611
Isotype	IgG
Reactivity	Human
Applications	WB;ELISA
Gene Name	KLHL25
Immunogen	Purified recombinant fragment of human KLHL25 expressed in E. Coli.
Specificity	KLHL25 Monoclonal Antibody detects endogenous levels of KLHL25 protein.
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1D8B8;Ectoderm-neural cortex protein 2;ENC 2;ENC2;ENC-2;KLHL25;Kelch-like protein 25
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm,Cul3-RING ubiquitin ligase complex,
Tissue Specificity	Brain,Testis,
Function	similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 6 Kelch repeats.,
Background	similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 6 Kelch repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using KLHL25 Monoclonal Antibody against KLHL25-hlgGfc transfected HEK293 cell.