



BM88 Polyclonal Antibody

Catalog No	YP-Ab-12690
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC;IF;ELISA
Gene Name	CEND1
Immunogen	The antiserum was produced against synthesized peptide derived from human CEND. AA range:11-60
Specificity	BM88 Polyclonal Antibody detects endogenous levels of BM88 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/5000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BM88 antigen;BM88;CEND1;Cell cycle exit and neuronal differentiation protein 1
Observed Band	23 kDa
Calculated Molecular Weight	149 aa, 15 kDa
Cell Pathway	Membrane ; Single-pass type IV membrane protein .
Tissue Specificity	Neuron specific.
Function	function:Involved in neuroblastoma cell differentiation.,similarity:Belongs to the CEND1 family.,subunit:Homodimer.,tissue specificity:Neuron specific.,
Background	The protein encoded by this gene is a neuron-specific protein. The similar protein in pig enhances neuroblastoma cell differentiation in vitro and may be involved in neuronal differentiation in vivo. Multiple pseudogenes have been reported for this gene. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images