




Cerebellin 2 Polyclonal Antibody


Catalog No	YP-Ab-12699
Isotype	IgG
Reactivity	Human;Mouse;Rat;Monkey
Applications	WB;ELISA
Gene Name	CBLN2
Immunogen	The antiserum was produced against synthesized peptide derived from human CBLN2. AA range:161-210
Specificity	Cerebellin 2 Polyclonal Antibody detects endogenous levels of Cerebellin 2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CBLN2;Cerebellin 2;cerebellin 2 precursor;Cerebellin-2
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Secreted .
Tissue Specificity	Brain,
Function	similarity:Contains 1 C1q domain.,subunit:May interact with CBLN1, CBLN3 and CBLN4.,
Background	similarity:Contains 1 C1q domain.,subunit:May interact with CBLN1, CBLN3 and CBLN4.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Cerebellin 2 Polyclonal Antibody



Western blot analysis of lysates from COS cells, using CBLN2 Antibody. The lane on the right is blocked with the synthesized peptide.