



Choactase Polyclonal Antibody

Catalog No	YP-Ab-12703
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	CHAT
Immunogen	The antiserum was produced against synthesized peptide derived from human Choactase. AA range:334-383
Specificity	Choactase Polyclonal Antibody detects endogenous levels of Choactase protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CHOACTase;Choline acetylase;choline acetyltransferase;Choline O-acetyltransferase;EC:2.3.1.6;CHAT
Observed Band	82,70kD
Calculated Molecular Weight	82,70kD
Cell Pathway	nucleus,cytoplasm,cytosol,presynapse,
Tissue Specificity	Brain,Lymphocyte,Placenta,Spinal cord,
Function	catalytic activity:Acetyl-CoA + choline = CoA + O-acetylcholine.,disease:Defects in CHAT are the cause of congenital myasthenic syndrome with episodic apnea (CMSEA) [MIM:254210]; formerly known as familial infantile myasthenia gravis 2 (FIMG2). CMSEA is an autosomal recessive congenital myasthenic syndrome. Patients have myasthenic symptoms since birth or early infancy, negative tests for anti-AChR antibodies, and abrupt episodic crises with increased weakness, bulbar paralysis, and apnea precipitated by undue exertion, fever, or excitement., function:Catalyzes the reversible synthesis of acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses.,online information:Choline acetyltransferase entry,similarity:Belongs to the carnitine/choline



acetyltransferase family.,

Background

This gene encodes an enzyme which catalyzes the biosynthesis of the neurotransmitter acetylcholine. This gene product is a characteristic feature of cholinergic neurons, and changes in these neurons may explain some of the symptoms of Alzheimer's disease. Polymorphisms in this gene have been associated with Alzheimer's disease and mild cognitive impairment. Mutations in this gene are associated with congenital myasthenic syndrome associated with episodic apnea. Multiple transcript variants encoding different isoforms have been found for this gene, and some of these variants have been shown to encode more than one isoform. [provided by RefSeq, May 2010],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Choactase Polyclonal Antibody diluted at 1:1000



Western blot analysis of lysate from HepG2 cells, using Choactase antibody.