




GRIN3 Polyclonal Antibody

Catalog No	YP-Ab-12735
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	GPRIN3
Immunogen	The antiserum was produced against synthesized peptide derived from human GPRIN3. AA range:131-180
Specificity	GRIN3 Polyclonal Antibody detects endogenous levels of GRIN3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPRIN3;KIAA2027;GRIN3;GPRIN family member 3;G protein-regulated inducer of neurite outgrowth 3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Brain,
Function	function:May be involved in neurite outgrowth.,
Background	function:May be involved in neurite outgrowth.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using GPRIN3 Antibody. The picture on the right is blocked with the synthesized peptide.

