




NTT5 Polyclonal Antibody

Catalog No	YP-Ab-12783
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	SLC6A16
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC6A16. AA range:233-282
Specificity	NTT5 Polyclonal Antibody detects endogenous levels of NTT5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NTT5;SLC6A16;Orphan sodium- and chloride-dependent neurotransmitter transporter NTT5;Solute carrier family 6 member 16
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Highly expressed in peripheral tissues, particularly in testis, pancreas, and prostate.
Function	similarity:Belongs to the sodium:neurotransmitter symporter (SNF) family.,tissue specificity:Highly expressed in peripheral tissues, particularly in testis, pancreas, and prostate.,
Background	SLC6A16 shows structural characteristics of an Na(+)- and Cl(-)-dependent neurotransmitter transporter, including 12 transmembrane (TM) domains, intracellular N and C termini, and large extracellular loops containing multiple N-glycosylation sites.[supplied by OMIM, Mar 2008],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemistry analysis of paraffin-embedded human brain, using SLC6A16 Antibody. The picture on the right is blocked with the synthesized peptide.