



SLC6A15 Polyclonal Antibody

Catalog No	YP-Ab-12807
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC
Gene Name	SLC6A15
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC6A15. AA range:32-81
Specificity	SLC6A15 Polyclonal Antibody detects endogenous levels of SLC6A15 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLC6A15;Sodium-dependent neutral amino acid transporter B(0)AT2;B0AT2;NTT73;SBAT1;Sodium- and chloride-dependent neurotransmitter transporter NTT73;Sodium-coupled branched-chain amino-acid transporter 1;Solute carrier family 6 member
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Almost exclusively expressed in the brain.
Function	function:Not yet known, orphan transporter.,similarity:Belongs to the sodium:neurotransmitter symporter (SNF) family.,
Background	This gene encodes a member of the solute carrier family 6 protein family which transports neutral amino acids. The encoded protein is thought to play a role in neuronal amino acid transport (PMID: 16185194) and may be associated with major depression (PMID: 21521612). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]




matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from Jurkat cells, using SLC6A15 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).