



FPRL2 Polyclonal Antibody

Catalog No	YP-Ab-13251
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IF;ELISA
Gene Name	FPR3
Immunogen	The antiserum was produced against synthesized peptide derived from human FPRL2. AA range:304-353
Specificity	FPRL2 Polyclonal Antibody detects endogenous levels of FPRL2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	fMet Leu Phe receptor;FMLP;fMLP receptor;formyl peptide receptor 1;FPR;FPRL2;FPR3;N-formyl peptide receptor 3;FPRH1;FMLP-related receptor II;FMLP-R-II;Formyl peptide receptor-like 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Detected in various tissues with highest expression in lung.
Function	function:Low affinity receptor for N-formyl-methionyl peptides, which are powerful neutrophils chemotactic factors. Binding of FMLP to the receptor causes activation of neutrophils. This response is mediated via a G-protein that activates a phosphatidylinositol-calcium second messenger system., similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	function:Low affinity receptor for N-formyl-methionyl peptides, which are powerful neutrophils chemotactic factors. Binding of FMLP to the receptor causes activation of neutrophils. This response is mediated via a G-protein that



activates a phosphatidylinositol-calcium second messenger system.,
similarity:Belongs to the G-protein coupled receptor 1 family.,



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunofluorescence analysis of LOVO cells, using FPRL2 Antibody. The picture on the right is blocked with the synthesized peptide.