



GPR116 Polyclonal Antibody

Catalog No	YP-Ab-13289
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC;IF;ELISA
Gene Name	GPR116
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR116. AA range:11-60
Specificity	GPR116 Polyclonal Antibody detects endogenous levels of GPR116 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADGRF5;Adhesion G protein-coupled receptor F5;G protein coupled receptor 116;G-protein coupled receptor 116;KIAA0758;GPR116;Probable G-protein coupled receptor 116
Observed Band	150kD
Calculated Molecular Weight	150kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain,Plasma,Testis,
Function	function:May have a role in the regulation of acid-base balance., PTM:Proteolytically cleaved into 2 highly conserved sites: one in the SEA domain and the other in the stalk domain region preceding the first transmembrane. The later 2 subunits, the extracellular subunit and the seven-transmembrane subunit, remain tightly associated and non-covalently linked., sequence caution:Contaminating sequence. Potential poly-A sequence., similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily., similarity:Contains 1 GPS domain.,similarity:Contains 1 SEA domain.,

similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,subunit:Exists as disulfide-linked dimers at the cell surface.,

Background

function:May have a role in the regulation of acid-base balance.,
PTM:Proteolytically cleaved into 2 highly conserved sites: one in the SEA domain and the other in the stalk domain region preceding the first transmembrane. The later 2 subunits, the extracellular subunit and the seven-transmembrane subunit, remain tightly associated and non-covalently linked.,
sequence caution:Contaminating sequence. Potential poly-A sequence.,
similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,
similarity:Contains 1 GPS domain.,similarity:Contains 1 SEA domain.,
similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,subunit:Exists as disulfide-linked dimers at the cell surface.,

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images




Western Blot analysis of various cells using GPR116 Polyclonal Antibody diluted at 1:500




Immunofluorescence analysis of HeLa cells, using GPR116 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from A549 cells, using GPR116 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using GPR116 antibody.



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).