



LDLRAD1 Polyclonal Antibody


Catalog No	YP-Ab-13390
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA;IHC
Gene Name	LDLRAD1
Immunogen	The antiserum was produced against synthesized peptide derived from human LDLRAD1. AA range:84-133
Specificity	LDLRAD1 Polyclonal Antibody detects endogenous levels of LDLRAD1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LDLRAD1;Low-density lipoprotein receptor class A domain-containing protein 1
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the LDLR family.,similarity:Contains 3 LDL-receptor class A domains.,
Background	similarity:Belongs to the LDLR family.,similarity:Contains 3 LDL-receptor class A domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Western Blot analysis of various cells using LDLRAD1 Polyclonal Antibody



Western blot analysis of lysates from Jurkat and HepG2 cells, using LDLRAD1 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).