



Cleaved-Caspase-4 p20 (Q81) Polyclonal Antibody

Catalog No	YP-Ab-00017
Isotype	IgG
Reactivity	Human
Applications	WB;IHC;IF;ELISA
Gene Name	
Immunogen	
Specificity	
Formulation	
Source	
Purification	
Dilution	
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cleaved-Caspase-4 p20 (Q81)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images