



NSD1 Monoclonal Antibody

Catalog No	YP-Ab-00089
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	NSD1
Immunogen	Purified recombinant human NSD1 protein fragments expressed in E.coli.
Specificity	NSD1 Monoclonal Antibody detects endogenous levels of NSD1 protein.
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Androgen receptor-associated protein of 267 kDa;ARA267;EC:2.1.1.357;H3-K36-HMTase;NSD1;Histone-lysine N-methyltransferase H3 lysine-36 and H4 lysine-20 specific;KMT3B;Histone-lysine N-methyltransferase;H3 lysine-36 and H4 lysine-20 specific;Androgen receptor coactivator 267 kDa protein;H4-K20-HMTase;Lysine N-methyltr
Observed Band	297kD
Calculated Molecular Weight	
Cell Pathway	Nucleus. Chromosome .
Tissue Specificity	Expressed in the fetal/adult brain, kidney, skeletal muscle, spleen, and the thymus, and faintly in the lung.
Function	catalytic activity:S-adenosyl-L-methionine + histone L-lysine = S-adenosyl-L-homocysteine + histone N(6)-methyl-L-lysine.,disease:A chromosomal aberration involving NSD1 is found in an adult form of myelodysplastic syndrome (MDS). Insertion of NUP98 into NSD1 generates a NUP98-NSD1 fusion product.,disease:A chromosomal aberration involving NSD1 is found in childhood acute myeloid leukemia. Translocation t(5;11)(q35;p15.5) with NUP98.,disease:Defects in NSD1 are a cause of Beckwith-Wiedemann syndrome (BWS) [MIM:130650]. BWS is a genetically heterogeneous disorder



characterized by anterior abdominal wall defects including exomphalos (omphalocele), pre- and postnatal overgrowth, and macroglossia. Additional less frequent complications include specific developmental defects and a predisposition to embryonal tumors. Disease: Defects in NSD1 are the cause of Sotos syndrome [MIM:117550]; also kn

Background

This gene encodes a protein containing a SET domain, 2 LXXLL motifs, 3 nuclear translocation signals (NLSs), 4 plant homeodomain (PHD) finger regions, and a proline-rich region. The encoded protein enhances androgen receptor (AR) transactivation, and this enhancement can be increased further in the presence of other androgen receptor associated coregulators. This protein may act as a nucleus-localized, basic transcriptional factor and also as a bifunctional transcriptional regulator. Mutations of this gene have been associated with Sotos syndrome and Weaver syndrome. One version of childhood acute myeloid leukemia is the result of a cryptic translocation with the breakpoints occurring within nuclear receptor-binding Su-var, enhancer of zeste, and trithorax domain protein 1 on chromosome 5 and nucleoporin, 98-kd on chromosome 11. Two transcript variants encoding distinct isofo


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis using NSD1 Monoclonal Antibody against HepG2 cell lysate .