



# CSRNP3 Polyclonal Antibody

|                             |                                                                                                                                                                                                                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                  | YP-Ab-00362                                                                                                                                                                                                                                                                       |
| Isotype                     | IgG                                                                                                                                                                                                                                                                               |
| Reactivity                  | Human;Mouse;Monkey                                                                                                                                                                                                                                                                |
| Applications                | WB;ELISA                                                                                                                                                                                                                                                                          |
| Gene Name                   | CSRNP3                                                                                                                                                                                                                                                                            |
| Immunogen                   | The antiserum was produced against synthesized peptide derived from human TAIP-2. AA range:171-220                                                                                                                                                                                |
| Specificity                 | CSRNP3 Polyclonal Antibody detects endogenous levels of CSRNP3 protein.                                                                                                                                                                                                           |
| Formulation                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                           |
| Source                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                                                                            |
| Purification                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                             |
| Dilution                    | Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.                                                                                                                                                                                               |
| Concentration               | 1 mg/ml                                                                                                                                                                                                                                                                           |
| Purity                      | ≥90%                                                                                                                                                                                                                                                                              |
| Storage Stability           | -20°C/1 year                                                                                                                                                                                                                                                                      |
| Synonyms                    | CSRNP3;Cysteine/serine-rich nuclear protein 3;FAM130A2;TAIP2;CSRNP-3; Protein FAM130A2;TGF-beta-induced apoptosis protein 2;TAIP-2                                                                                                                                                |
| Observed Band               | 64kD                                                                                                                                                                                                                                                                              |
| Calculated Molecular Weight | 64kD                                                                                                                                                                                                                                                                              |
| Cell Pathway                | Nucleus .                                                                                                                                                                                                                                                                         |
| Tissue Specificity          | Embryo,Thymus,                                                                                                                                                                                                                                                                    |
| Function                    | function:Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis.,sequence caution:Unusual initiator. The initiator methionine is coded by a non-canonical CTG leucine codon.,similarity:Belongs to the AXUD1 family., |
| Background                  | function:Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis.,sequence caution:Unusual initiator. The initiator methionine is coded by a non-canonical CTG leucine codon.,similarity:Belongs to the AXUD1 family., |


**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images



Western blot analysis of the lysates from HepG2 cells using TAIP-2 antibody.